

# BDC-600

## Dry Biochemical analysis meter

Take your lab's performance to the next level



### Blood Lipid Item:

- » Total Cholesterol (TC)
- » High Density Lipoprotein (HDL)
- » Triglyceride (TG)

### Renal Function Item:

- » Uric Acid (UA)
- » Creatinine (CR)
- » Urea (UR)

## Advantage

### Multi-device function

One meter with 6 parameters

### Wide Range

Maximum measurement range in same function

POCT meter

### Quality Assurance

Excellent quality and advanced technology

### Accurate

Coincide with the result from Biochemical Analyzer

### Portable

Light and small enough to be easily carried or moved

### Rechargeable

Built-in 1200mAh Li-ion battery and support

Long-Endurance



## Specification

Methodology	Reflectance Photometer
Specimen	Whole Blood - Serum - plasma
Power Source	1200mAh built-in lithium battery
Battery Life	After cycle charge 300 times, electric quantity reduce 30%
Units of Measurement	mg/dL, mmol/L
Memory	500 records
Automatic Shut Off	5 minutes after last use
Meter Size	135mm X 66mm X 19mm
Weight	90g
Meter Storage Conditions	0-55 °c – 90%RH
System Operating Conditions	10-35 °c – 90%RH: height – seal level 2000m
Test strip Storage Conditions	2-30 °c – 90%RH
Warranty Period	5 years

Inzek B.V.

Vissenstraat 32 | 7324 AL Apeldoorn | The Netherlands | info@biozek.com | www.biozek.com